



Med Instrument Trading llc. is a dynamic and innovative medical instrument trading company dedicated to sourcing, distributing, and supplying high-quality medical equipment and instruments to healthcare facilities worldwide. With a commitment to excellence, reliability, and customer satisfaction, we have established ourselves as a reliable and reputable partner in the medical industry.

Med Instrument Medical Equipment Trading LLC

Duabi Investment Park 1, Falcon House
Second Floor, Czar Business Center
Section 204 and 203

Tel: 00971 55 10 2525 1
Email: info@medinstrument.net
Web: www.medinstrument.net



SUPER CUT SCISSORS


Manufactured from the Highest Quality Stainless Steel

The Supercut Advantage

All Supercut Scissors are expertly manufactured from the highest quality stainless steel. The result is a significantly sharper instrument with super fine precision capabilities, making them ideal for a variety of surgical procedures.

1. What are Supercut Scissors ?

Supercut Scissors offer unrivalled quality and cutting ability and are "sharper" than conventional scissors. The cutting blade is "razor-edged".



2. Why Supercut Scissors are Better than the Ordinary Scissors?

A specially designed razor-sharp upper blade edge that cuts effortlessly through tissue. Combined with a lower blade having micro serrations to hold and prevent tissue slippage make these scissors ideal for gripping tissue and cutting interproximally.



3. What are their Other Advantages?

Supercut Scissors with razor edge + micro serrations ensure extremely controlled cutting in difficult areas during surgery. They cut their way forward more accurately.

These scissors supersede because of their ability to cut with extreme precision right down to the tip of the blades.

Literally speaking, the cutting edges obey the hands and mind of the surgeon most precisely. The cutting margins/ edges are honed to eliminate any traumatization of the tissues during cutting action, by applying special custom-made technology.

For easy identification all Hipp and Armo Supercut Scissors have black ring handles.



LISTER BANDAGE SCISSORS

Lister Bandage Scissors are used for trimming dressings and removing bandages as well as cutting tougher material such as plaster and fabric.

14.5cm	M-S1-01
18cm	M-S1-02
20cm	M-S1-03



UNIVERSAL SCISSORS
BLACK



UNIVERSAL SCISSORS
ORANGE



UNIVERSAL SCISSORS
PINK



UNIVERSAL SCISSORS
GREEN



UNIVERSAL SCISSORS
BLUE



UNIVERSAL SCISSORS
PURPLE

UNIVERSAL SCISSORS

Universal/Trauma scissors are used by paramedics and other emergency medical personnel to cut clothing. They have a large finger ring for easy handling and angled serrated blade for precise cutting. They are autoclavable and have coloured handles and plastic tip.

LARGE	19cm
Black	M-US2-01
Green	M-US2-02
Red	M-US2-03
Pink	M-US2-04

SMALL	16cm
Black	M-US2-05
Green	M-US2-06
Red	M-US2-07
Blue	M-US2-08
Orange	M-US2-09
Pink	M-US2-10
Purple	M-US2-11
Yellow	M-US2-12



WARD SCISSORS


Ward scissors are general purpose scissors commonly used for cutting bandages and gauze, available with or without black rubber coated handle.

20cm	M-S3-01
20cm	M-S3-02

**IRIS SCISSORS
TUNGSTEN CARBIDE**

Iris scissors are surgical scissors with short sharp tips. Ideal for precision tasks and ophthalmic applications.



Straight	11cm	M-TC8-01
Curved	11cm	M-TC8-02



CUTICLE SCISSORS

Cuticle scissors are finely pointed and thin with a small blade in order to trim tiny hangnails or frayed tissue around the nail.


Straight	9cm	M-S4-01
Curved	9cm	M-S4-02



**IRIS SUPERCUT
SCISSORS**

Iris scissors are surgical scissors with short sharp tips. Ideal for precision tasks and ophthalmic applications.

Straight	11cm	M-SC9-01
Curved	11cm	M-SC9-02



NAIL SCISSORS



Nail scissors are small scissors with a slender shaft and strong blades that are used for shaping and trimming fingernails.

Straight	9cm	M-S5-01
Curved	9cm	M-S5-02

KILNER SCISSORS

Kilner scissors are used for the dissecting and cutting of delicate tissue.


Straight	13cm	M-S10-01
Curved	13cm	M-S10-02



CASTROVIEJO SCISSORS

Castroviejo scissors are double spring ophthalmic scissors commonly used during iridectomy and corneal surgery.


Straight	20cm	M-S6-01
Curved	20cm	M-S6-02



**KILNER SUPERCUT
SCISSORS**

Kilner scissors are used for the dissecting and cutting of delicate tissue.

Curved	13cm	M-SC11-01
--------	------	-----------



IRIS SCISSORS


Iris scissors are surgical scissors with short sharp tips. Ideal for precision tasks and ophthalmic applications.



Straight	20cm	M-S7-01
Curved	20cm	M-S7-02

STRABISMUS SCISSORS

Strabismus scissors are blunt tipped scissors used for ophthalmic dissection.

Straight	11cm	M-S12-01
Curved	11cm	M-S12-02





STRAIGHT CURVED

KELLY SCISSORS

Kelly scissors are sharp tipped and are used to grasp and dissect fascia and soft tissue.

Straight	16cm	M-S13-01
Curved	16cm	M-S13-02

FERGUSSON SCISSORS

Kelly scissors are sharp tipped and are used to grasp and dissect fascia and soft tissue.

18cm	M-S17-01
21cm	M-S17-02



CURVED

STEVENS SCISSORS

Stevens Tenotomy Scissors are used for delicate dissection and cutting, usually in ophthalmic, neurological and plastic surgery procedures.

Straight	11cm	M-S14-01
Curved	11cm	M-S14-02



STRAIGHT CURVED



BRAUN STADLER SCISSORS

Braun-Stadler Episiotomy Scissors are used to perform episiotomies during labour and delivery.

14.5cm	M-S18-01
22cm	M-S18-02



STRAIGHT CURVED

STEVENS SUPERCUT SCISSORS

Stevens Tenotomy Scissors are used for delicate dissection and cutting, usually in ophthalmic, neurological, and plastic surgery procedures.

Straight	11.5cm	M-SC15-01
Curved	11.5cm	M-SC15-02

UMBILICAL SCISSORS

Umbilical scissors used for cutting the umbilical cord.


10.5cm	M-S19-01
12.5cm	M-S19-02



O'BRIEN SCISSORS

O'Brien Scissors are primarily used for suture removal

9cm	M-S16-01
11.5cm	M-S16-02



STRAIGHT CURVED


SIMS UTERINE SCISSORS

Sims (Uterine) scissors have long blades, primarily used for deep tissue dissection.

Straight	20cm	M-S20-01
Curved	20cm	M-S20-02
	23cm	M-S20-03

MAYO SCISSORS


Mayo Scissors have blunt ends and are generally used for cutting sutures & ligatures or trimming the skin around a wound.



Straight	14.5cm	M-S21-01	20cm	M-S21-04
	16cm	M-S21-02	22cm	M-S21-05
	18cm	M-S21-03		
Curved	14.5cm	M-S21-06	20cm	M-S21-09
	16cm	M-S21-07	22cm	M-S21-10
	18cm	M-S21-08		

METZENBAUM SCISSORS


Metzenbaum scissors are surgical scissors designed for cutting delicate tissue and blunt dissection. These scissors feature a more delicate profile than the standard metzenbaum scissor.



Straight	12cm	M-S25-01	18cm	M-S25-04
	14cm	M-S25-02	20cm	M-S25-05
	16cm	M-S25-03	22cm	M-S25-06
Curved	12cm	M-S25-07	20cm	M-S25-11
	14cm	M-S25-08	22cm	M-S25-12
	16cm	M-S25-09	28cm	M-S23-13
	18cm	M-S25-10		

MAYO SCISSORS TUNGSTEN CARBIDE



Mayo Scissors have blunt ends and are generally used for cutting sutures and ligatures or trimming the skin around a wound.



Straight	14.5cm	M-TC22-01
	18cm	M-TC22-02
	22cm	M-TC22-03
Curved	14.5cm	M-TC22-04
	18cm	M-TC22-05
	22cm	M-TC22-06

METZENBAUM SCISSORS FINE


Metzenbaum scissors are surgical scissors designed for cutting delicate tissue and blunt dissection. They are finer and more delicate than Mayo scissors.



Fine Tip Straight	22cm	M-S26-01
Fine Tip Curved	22cm	M-S26-02

MAYO SUPERCUT SCISSORS


Mayo Scissors have blunt ends and are generally used for cutting sutures and ligatures or trimming the skin around a wound.



Straight	14cm	M-SC23-01
	17cm	M-SC23-02
	20cm	M-SC23-03
Curved	14.5cm	M-SC23-04
	17cm	M-SC23-05

METZENBAUM SCISSORS TUNGSTEN CARBIDE



Tungsten carbide metzenbaum scissors are surgical scissors designed for cutting delicate tissue and blunt dissection. They are finer and more delicate than mayo scissors. Having the blades infused with tungsten carbide allows for greater durability of the scissor.



Straight	14.5cm	M-TC27-01	Curved	14.5cm	M-TC27-05
	18cm	M-TC27-02		18cm	M-TC27-06
	20cm	M-TC27-03		20cm	M-TC27-07
	22cm	M-TC27-04		22cm	M-TC27-08

MAYO HARRINGTON SCISSORS


Mayo Harrington Dissecting Scissors have flat sharp blades with blunt tips and are available in straight or curved in a 22cm length.



Straight	22cm	M-S24-01
Curved	22cm	M-S24-02

METZENBAUM SUPERCUT SCISSORS

Supercut metzenbaum scissors are surgical scissors designed for cutting delicate tissue and blunt dissection. They are finer and more delicate than Mayo scissors.





Straight	9cm
Curved	9cm

Straight	14cm	M-SC28-01	Curved	14cm	M-SC28-03		
	Sharp Blunt	18cm		M-SC28-02	Sharp Blunt	18cm	M-SC28-04
Straight	11cm	M-SC28-05	Curved	11cm	M-SC28-08		
	Blunt/Blunt	14cm		M-SC28-06	Blunt/Blunt	14cm	M-SC28-09
	18cm	M-SC28-07		18cm	M-SC28-10		

SURGICAL SCISSORS

Surgical or operating scissors are used for general cutting and dissecting of tissue, bandages and gauze.



STRAIGHT

Straight	14cm	M-S29-01	20cm	M-S29-04
Sharp/Sharp	16cm	M-S29-02	22cm	M-S29-05
	18cm	M-S29-03	24cm	M-S29-06

Straight	14cm	M-S29-07	20cm	M-S29-10
Sharp/Blunt	16cm	M-S29-08	22cm	M-S29-11
	18cm	M-S29-09		

Straight	14cm	M-S29-12	20cm	M-S29-15
Blunt/Blunt	16cm	M-S29-13	22cm	M-S29-16
	18cm	M-S29-14		


Straight	14cm	M-S29-17		
Sharp/Probe	16cm	M-S29-18		

CURVED

Curved	14cm	M-S29-19	20cm	M-S29-22
Sharp/Sharp	16cm	M-S29-20	22cm	M-S29-23
	18cm	M-S29-21	24cm	M-S29-24

Curved	14cm	M-S29-25	20cm	M-S29-28
Sharp/Blunt	16cm	M-S29-26	22cm	M-S29-29
	18cm	M-S29-27		

Curved	14cm	M-S29-30	20cm	M-S29-33
Blunt/Blunt	16cm	M-S29-31	22cm	M-S29-34
	18cm	M-S29-32		



SPENCER SUTURE SCISSORS

Spencer Scissors are primarily used for suture removal. These scissors have a hook shaped tip on one blade that slides under sutures to slightly lift them before cutting.

9cm	M-S31-01
12cm	M-S31-02
14cm	M-S31-03



GOLDMANN FOX SUPERCUT SCISSORS

Goldmann Fox Scissors are used in dental procedures to trim gingival tissue. They are also used in plastic surgery or hemorrhoidectomy. They have curved, bevelled blades, one of which is serrated, that taper into sharp, fine tips.

Curved	13cm	M-SC32-01
--------	------	-----------



SURGICAL SCISSORS TUNGSTEN CARBIDE

Tungsten carbide surgical or operating scissors are used for general cutting and dissecting of tissue, bandages and gauze. They are available with blunt or sharp blades or a combination of both.

STRAIGHT

Straight	14cm	M-TC30-01
Sharp/Sharp	16cm	M-TC30-02
	18cm	M-TC30-03

Straight	14cm	M-TC30-04
Sharp/Blunt	16cm	M-TC30-05
	18cm	M-TC30-06

Straight	14cm	M-TC30-07
Blunt/Blunt	16cm	M-TC30-08
	18cm	M-TC30-09

CURVED

Curved	14cm	M-TC30-10
Sharp/Sharp	16cm	M-TC30-11
	18cm	M-TC30-12

Curved	14cm	M-TC30-13
Sharp/Blunt	16cm	M-TC30-14
	18cm	M-TC30-15

Curved	14cm	M-TC30-16
Blunt/Blunt	16cm	M-TC30-17
	18cm	M-TC30-18



REYNOLDS SUPERCUT SCISSORS

Reynolds scissors are characterised by having narrow, tenotomy style tips, generally used for cutting tissue as well as dissection.

Curved	15cm	M-SC33-01
--------	------	-----------

WAGNER SUPERCUT SCISSORS

Wagner Scissors are used in procedures requiring delicate or precise cutting, particularly plastic surgery. Available in 12cm length with straight or curved Supercut blades

Straight	12cm	M-SC34-01
Curved	12cm	M-SC34-02

